

Release Notes

ACH MEDICAL Digital Library
(ver1.1)



Development Team



5 Page



-



2026.01.09

Table of Contents



1. Release Notes	3
2. Reference Image.....	4
2-1. EOXCAD	
2-2. 3SHAPE	
3. Revision Summary	5

1. Release Notes

2026.01.09

Applying Client Requests (25.11.23)

- Verification and Correction of Multi-Unit Library TiBase Code Errors (ex. 4100 → 0410)



- Newly Added Multi-Unit Screw-Retained Type Library

- The initial library was produced only for the Cylinder Type (Multi-TiBase)

: (Method of fabricating crown/bridge after fixing Multi-TiBase on Multi Abutment with cement or resin)


<NEW>

- Newly created library for Screw-Retained Type (Multi-Unit direct screw-retained prosthesis)


: Fabrication method where the prosthesis (Crown) itself is directly screwed onto the Multi Abutment) → No cylinder used

2. Reference Image

2-1. EXOCAD



2-2. 3Shape



3. Revision Summary

Version	Date	Summary of Changes	Author
1.1	2026.01.09	Added Multi-Unit Screw-Retained Type Only the Multi-Unit library has been updated The corresponding library file: ① ACH Medical EXOCAD integrated Library_ver1.1(26.01.09) ② ACH Medical 3Shape integrated Library_ver1.1(26.01.09)	Development Team
1.0	2025.09.30	- Re-validation of previous library data - Added IB MINI Green Digital Lineup Library - Improved and modified S/W file names for better user recognition and selection convenience - Integrated construction for library internalization	Development Team